

OXID MICROMANIPULATOR

ADAPTER CHART

12/19/2005

PART #	DESCRIPTION
PLADR0068	OXID DOVETAIL TO ZEISS 65mm DIA CLAMP
PLADR0072	OXID DOVETAIL TO WILD M62 THREAD (FINE PITCH)
PLADR0074	OXID DOVETAIL TO CODMAN M55 THREAD
PLADR0199	OXID DOVETAIL TO ZEISS M47 THREAD
PLADR0315	OXID DOVETAIL TO LEICA M75 THREAD
PLADR0321	OXID DOVETAIL TO WILD M62 THREAD (COARSE PITCH)
PLADR0338	OXID DOVETAIL TO NIKON M48 THREAD
PLADR0384	OXID DOVETAIL TO ZEISS DOVETAIL